



DOUBLE PROBE ULTRASOUND - D5CE

AFFORDABLE
VISIBLE - FAST - EASY



YOU CAN USE IT ANYTIME & ANYWHERE

UNIQUE MEDI TRADE

Begum Rokeya Sharani, Shre-E-Bangla Nagar, Dhaka

E-mail: contact@uniquemeditrade.com

Telephone: +880 2 55025199

Phone/WhatsApp: +880 1717 81 13 12

Website: www.uniquemeditrade.com



APPLICATION

Application Sites:

- OB/GYN
- Urology

Application Directions:

Follicular Monitoring, Prostate/Pelvics, Uterus, Bladder/Urinary Function.



FEATURES

- Composed of two probes, a scanner can complete surface and vaginal examination as your choice.
- Small and light, suitable for bedside operation and monitoring.
- Wireless, by using the disposable protecting cover, it can easily solve the sterilization issue of the probe.



SPECIFICATION

- Scan Angle/Length: Convex 60, Linear 40mm
- Frequency: Convex probe 3.5/5.0 MHz, Linear probe 6.0/8 MHz
- Scan Mode: B, B/M, Color, PW
- Display Depth: Convex 90~305mm, Transvaginal 40~100mm
- Power: Built-in lithium battery, wireless charge 2h, work 2.5h
- Working System: iOS, Android
- WiFi Type: 802.11n / 2.4G / 5G
- Dimension: 270x60x20mm
- Weight: 300g



WATERPROOF

Reduce medical accidents caused by cross infection, protect doctors and patients, and reduce hospital related expenses.



WIRELESS CHARGE

Easy and connect with a wireless power bank.



WiFi CONNECT

WiFi connect with your Android phone, iPhone or iPad.



EASY TO USE

Simple operation, short learning curve for people without ultrasound experience.



HIGH QUALITY IMAGES

See The Difference Great Image Quality Makes

