



# WIRELESS PROBE ULTRASOUND - P5L

**AFFORDABLE**  
VISIBLE - FAST - EASY



YOU CAN USE IT ANYTIME & ANYWHERE

## UNIQUE MEDI TRADE

Begum Rokeya Sharani, Shre-E-Bangla Nagar, Dhaka  
 E-mail: [contact@uniquemeditrade.com](mailto:contact@uniquemeditrade.com)  
 Telephone: +880 2 55025199  
 Phone/WhatsApp: +880 1717 81 13 12  
 Website: [www.uniquemeditrade.com](http://www.uniquemeditrade.com)

### APPLICATION

- Application Sites: Small Parts, MSK, Nerve, Vascular access  
 Application Directions:
- Medical Students training in college: Visual understanding of human body structure or ultrasound image study.
  - Preliminary examination: Home care, family planning.
  - Visualization tools: Invasive intervention guide, surgical and therapy guidance.
  - ODM/OEM: Integration into other multi-functional medical device.

### FEATURES

- Super Cost-effective: Students can afford it, just like buying a cell phone.
- Small, light, built in battery: You can easily use it anywhere.
- Remote training and guidance: Based on intelligent terminal, remote training and guidance can be achieved.

### SPECIFICATION

- Scan Angle/Length: 40mm
- Frequency: 7.5/10 MHz
- Scan Mode: B, B/M, Color, PW
- Display Depth: 20~100mm
- Power: Built-in lithium battery, USB charge 1.5h, work 1.5h
- Working System: iOS, Android
- WiFi Type: 802.11n / 2.4G / 5G
- Dimension: 115x60x20mm
- Weight: 100g



### USB CHARGE

Just simply connect to your laptop or USB charger.



### EASY TO USE

Simple operation, short learning curve for people without ultrasound experience.



### WiFi & USB CONNECT

- Both WiFi and USB connect as you can choose.
- Use with your Android phone, Laptop, iPhone or iPad.



### HIGH QUALITY IMAGES

See The Difference Great Image Quality Makes

