



Handheld Ultrasound  
Scanner Manufacturer



# WIRELESS ULTRASOUND

## AFFORDABLE VISIBLE - FAST - EASY

**GUANGZHOU SONOHEALTH MEDICAL  
TECHNOLOGIES CO., LTD.**

Address: Room 601, Building 3, No.7 Ruitai Road, Huangpu District, Guangzhou, China

E-mail: [hello@sono-health.com](mailto:hello@sono-health.com)

Telephone: +86-02-87831007

Phone/WhatsApp: +86-15013156300

Website: [www.sono-health.com](http://www.sono-health.com)



YOU CAN USE IT ANYTIME & ANYWHERE



## APPLICATION

- Application Sites
- OB/GYN
- Application Directions
- Follicular monitoring



## FEATURES

- Low price, individuals can afford to achieve home monitoring, avoid frequent hospital visits.
- Small and light, suitable for bedside and home monitoring.
- Wireless, by using the disposable protecting cover, it can easily solve the sterilization issue of the probe.



## SPECIFICATION

- Scan Angle/Length: 149
- Frequency: 6/8 MHz
- Scan Mode: B, B/M, Color, PW
- Display Depth: 30~115mm
- Power: Built-in lithium battery, USB charge 1.5h, work 1.5h
- Working System: iOS, Android
- WiFi Type: 802.11n / 2.4G / 5G
- Dimension: 336x35x28mm
- Weight: 170g



## WIRELESS & USB CHARGE

- Both wireless charge and USB charge you can choose.
- Easy and connect with a wireless power bank or USB charger.

## EASY TO USE

Simple operation, short learning curve for people without ultrasound experience.



## WiFi & USB CONNECT

- Both WiFi and USB connect as you can choose.
- Use with your Android phone, Laptop, iPhone or iPad.

## HIGH QUALITY IMAGES

See The Difference Great Image Quality Makes

